

## In-situ Straight Leg Sling

### S17

- Designed to provide full body support on transfer
- Easy to apply with low-friction fabric to leg areas for ease of application
- Soft binding as standard
- Fixed colour coded loops or Silva-Safety Slot® tapes
- Once the patient has been transferred, this sling can be left in-situ, subject to an appropriate risk assessment
- Wash Temp 71°C/159°F
- Tumble dry low heat (max 50°C)



nb: Sling shown in image is with fixed Silva-Safety Slots® tapes

### Product Codes & Specifications:



	Silva-Superfine PLUS®	Silva-Superfine® with COOLMAX®	Silva-Superfine® with COOLMAX® Double Layer	A	B	C	D
XS	-	S17CXS	-	76 - 86cm	44 - 50cm	76cm	55cm
S	S17FSM	S17CS	-	86 - 96cm	50 - 57cm	80cm	58cm
M	S17FMED	S17CM	S17SFCDLMED	96 - 108cm	57 - 64cm	90cm	66cm
L	S17FLRG	-	S17SFCDLLRG	108 - 130cm	64 - 69cm	102cm	74cm
XL	S17FXLRG	-	S17SFCDLXLRG	130 - 144cm	69 - 74cm	120cm	78cm

