

## **In-situ Deluxe Leg Sling**

- Designed to provide full body support on transfer
- · Shaped leg sections provide full femur & hip support
- Easy to apply with low-friction fabric to leg areas for ease of application
- Soft binding as standard
- Fixed colour coded loops or Silva-Safety Slot® tapes
- Once the patient has been transferred, this sling can be left in-situ, subject to an appropriate risk assessment
- Wash Temp 71°C/159°F
- Tumble dry low heat (max 50°C)


## **Product Codes & Specifications:**










	Silva-Superfine PLUS®	Silva-Superfine® with COOLMAX®	Silva-Superfine® with COOLMAX® Double Layer	A	В	C	D
xxs	S18FXXS	S18CXXS	-	56 - 66cm	38 - 44cm	58cm	49cm
xs	S18FXS	S18CXS	-	76 - 86cm	44 - 50cm	68cm	54cm
s	S18FSM	S18CS	-	86 – 96cm	50 – 57cm	79cm	58cm
м	S18FMED	S18CM	S18CDL-M	96 – 108cm	57 – 64cm	90cm	65cm
L	S18FLRG	-	S18CDL-L	108 – 130cm	64 – 69cm	98cm	72cm
XL	S18FXLRG	-	S18CDL-XL	130 – 144cm	69 – 74cm	114cm	78cm
XXL	S18FXXL	-	S18CDL-XXL	144 - 160cm	69 – 74cm	122cm	82cm









Silvalea Ltd a division of **3 savaria**. Copyright 2022 **2** All rights reserved

## **S18**

