

Sit to Stand Sling

405

- The Sit to Stand harness is available with a front or reverse fastening buckle
- Ideal for transferring to a standing frame
- Used to transfer in a 90 degree position
- Useful to reposition patients legs once seated in wheelchair
- Can be used on an overhead hoist or a mobile hoist
- Wash Temp 85°C/185°F
- Tumble dry low heat (max 50°C)



Product Codes & Specifications:



100% Polyester

		A	B	C
XXS	405XXS	77 - 87cm	61cm	15cm
XS	405XS	87 - 107cm	75cm	17cm
S	405S	97 - 117cm	91cm	20cm
M	405M	112 - 132cm	99cm	25cm
L	405L	125 - 145cm	109cm	28cm
XL	405XL	138 - 158cm	127cm	31cm

