

## Silva Jolly Jumper

- Good Supportive harness
- Only two hoisting straps attached to the top of the head support allowing more freedom for the patient
- This product is supplied with adjustable sliders as standard
- Provides support to improve standing tolerance and independent transfers
- Used for rehabilitation purposes, therapy and inclusion in school & family activities
- Can be used on an overhead hoist or a mobile hoist
- Manufactured in quilted polyester
- Wash Temp 85°C/185°F. Tumble dry low heat (max 50°C)

## SILJJ



### Product Codes & Specifications:



| 100% Polyester |           | A         | B         | C    | D    |
|----------------|-----------|-----------|-----------|------|------|
| XXXS           | SILJJXXXS | 38-58cm   | 50 - 55cm | 24cm | 47cm |
| XXS            | SILJJXXS  | 59-77cm   | 56 - 62cm | 27cm | 77cm |
| XS             | SILJJXS   | 77-87cm   | 60 - 66cm | 33cm | 82cm |
| S              | SILJJ S   | 87-107cm  | 66 - 74cm | 38cm | 87cm |
| M              | SILJJM    | 97-117cm  | 74 - 80cm | 43cm | 91cm |
| L              | SILJJL    | 112-132cm | 80 - 88cm | 45cm | 91m  |

