


## Back and Leg Raiser

### 201

- Designed to assist carer with repositioning a patient
- Helps to reduce handling
- Can be used across shoulders to gain access to insert a sling
- Can be used to assist a patient into a sitting position when in bed or chair
- Can be used either on one leg or both together to raise or reposition and to assist with any wound care treatment
- Option of both inner and outer positioning handles
- This is a manual device and current moving and handling guidelines should be followed
- Manufactured in a quilted, polyester fabric with an inner soft lining for additional comfort
- Wash temp 71°C/159°F
- Tumble dry low heat (max 50°C)



### Product Codes & Specifications:

|                 |   |          |          |
|-----------------|---|----------|----------|
|                 |  |          |          |
|                 | <b>Polyester</b>  | <b>A</b> | <b>B</b> |
| <b>ONE SIZE</b> | <b>201</b>  | 58cm     | 22cm     |

