

In-situ Deluxe Leg Sling

S18

- Designed to provide full body support on transfer
- Shaped leg sections provide full femur & hip support
- Easy to apply with low-friction fabric to leg areas for ease of application
- Soft binding as standard
- Fixed colour coded loops or Silva-Safety Slot® tapes
- Once the patient has been transferred, this sling can be left in-situ, subject to an appropriate risk assessment
- Wash Temp 71°C/159°F
- Tumble dry low heat (max 50°C)



nb: Sling shown in image is with fixed Silva-Safety Slots® tapes

Product Codes & Specifications:



	Silva-Superfine PLUS®	Silva-Superfine® with COOLMAX®	Silva-Superfine® with COOLMAX® Double Layer	A	B	C	D
XXS	S18FXXS	S18CXXS	-	56 - 66cm	50 - 55cm	58cm	49cm
XS	S18FXS	S18CXs	-	76 - 86cm	60 - 66cm	68cm	54cm
S	S18FSM	S18CS	-	86 - 96cm	50 - 57cm	79cm	58cm
M	S18FMED	S18CM	S18CDL-M	96 - 108cm	57 - 64cm	90cm	65cm
L	S18FLRG	-	S18CDL-L	108 - 130cm	64 - 69cm	98cm	72cm
XL	S18FXLRG	-	S18CDL-XL	130 - 144cm	69 - 74cm	114cm	78cm
XXL	S18FXXL	-	S18CDL-XXL	144 - 160cm	96 - 104cm	122cm	82cm

