

## Trunk & Pelvic Walking Harness

## NWH023

- Provides significant support, particularly around the trunk and pelvic areas due to its elasticsation
- Helps to improve standing tolerance, independent transfers and activities using natural patterns of movement
- Successfully used with children with athetoid cerebral palsy, acquired brain injury and those with ataxic movement
- Allows forward and backward movement for balance activities
- Suitable for rehabilitation purposes, therapy, working in prone positions and inclusion in school and family activities
- Suitable for hydrotherapy purposes
- Needs to be applied on a bed or changing table
- Comfort covers can be added to the straps



### Product Codes & Specifications:



	Neoprene	A	B	C	D
XXS	NWH023XXSM	43 - 59cm	34 - 38cm	78cm	21cm
XS	NWH023XSM	59 - 77cm	38 - 44cm	88cm	24cm
XS	NWH023SM	77 - 87cm	44 - 50cm	98cm	30cm
M	NWH023MED	87 - 107cm	50 - 57cm	108cm	35cm
L	NWH023LRG	97 - 117cm	57 - 64cm	120cm	42cm

