# **Tape Fixing Options**




### **Universal Clip**

Codes: DA or CP

Our Universal D-anserclip® fits all the major hoist models within the UK and European markets, allowing greater interchangeability between varying clip hoists.

The safety "Click-on Click-off" action ensures a physical click is heard on the application of the spreader bar. All moving parts are internal to the clip ensuring greater durability during the lifetime of the patient sling.

#### **Compatible for Stud Clip Hoists only.**





### Silva-Safety Slots®

Codes: SSS


This unique, patented solution replaces the standard loops on patient transfer slings with a strong, durable, slotted webbing. Vibrant pick- up colours ensure the accuracy of application by caregivers.

The elimination of loops helps reduce risks associated with loops being caught in machinery (hoists or chairs) or creating a trip hazard when a sling is placed onto the floor.

The streamlined Silva-Safety Slots® ensure less bulk for those using in-situ slings where space is limited in wheelchairs and other seating options.

#### Compatible for Coat Hanger (Loop) Hoists only.





### **Codes: LP**


## **Colour Coded Loops**

Popular loop fixing option for standard patient transfer slings. Multiple colour coded internal loops offer maximum flexibility when positioning patients

This ensures optimum accuracy of application for caregivers as the coloured positions can be recorded as part of a care plan for consistancy.

#### Compatible for Coat Hanger (Loop) Hoists only.






For further advice, please refer to our full 'Sling & Hoist Compatibility Statement' found **HERE** or at **silvalea.com** on our download page.









Email: customerservices@silvalea.com

Website: silvalea.com