


Repositioning Sheet

600

- Used for repositioning and lateral transfers
- Used to help with effective pressure management, to reduce and prevent pressure sores
- Constructed with reinforced tapes from the side of the sling
- Can be used as part of the bed make-up when used in combination with a draw sheet
- Generous size making it suitable for all body dimensions
- Wash Temp 85°C/185°F
- Tumble dry low heat (max 50°C)



Product Codes & Specifications:

				
				
	100% Close Weave Polyester - Heavy Duty	100% Polyester Mesh	A	B
One Size	600A	600B	132cm	247cm

