

## In-situ High Easy All in One

## DA-HYA2

- Designed to provide full body & head support on transfer
- Apply on the bed or in a supine position
- Suitable for amputee scenarios
- Soft binding as standard
- Universal clip fixing
- Once the patient has been transferred, this sling can be left in-situ, subject to an appropriate risk assessment
- Wash Temp 71°C/159°F
- Tumble dry low heat (max 50°C)



### Product Codes & Specifications:



Silva-Superfine PLUS®		A	B	C	D
XXS	DA-HYA2-F-XXS	66 - 76cm	56 - 62cm	66cm	83cm
XS	DA-HYA2-F-XS	76 - 86cm	60 - 66cm	68cm	88cm
S	DA-HYA2-F-S	86 - 96cm	66 - 74cm	76cm	98cm
M	DA-HYA2-F-M	96 - 108cm	74 - 80cm	86cm	107cm
L	DA-HYA2-F-L	108 - 130cm	80 - 88cm	98cm	112cm
XL	DA-HYA2-F-XL	130 - 144cm	88 - 100cm	102cm	118cm
XXL	DA-HYA2-F-XXL	144 - 160cm	96 - 104cm	120cm	124cm

