

In-situ Deluxe Leg Sling HS

S16

- Designed to provide full body support on transfer
- Adjustable integral head support with removable stiffeners
- Shaped leg sections provide full femur & hip support
- Easy to apply with low-friction fabric to leg areas for ease of application
- Soft binding as standard
- Fixed colour coded loops
- Once the patient has been transferred, this sling can be left in-situ, subject to an appropriate risk assessment
- Wash Temp 71°C/159°F
- Tumble dry low heat (max 50°C)



nb: Sling shown in image is an In-Situ Deluxe Leg Sling HS with fixed Silva-Safety Slots*

Product Codes & Specifications:



| | Silva-Superfine PLUS® | Silva-Superfine® with COOLMAX® | Silva-Superfine® with COOLMAX® Double Layer | A | B | C | D |
|------|-----------------------|--------------------------------|---|-------------|------------|-------|-------|
| XXXS | S16FXXXSM | S16SFCXXXSM | S16CDLXXXS | 56 - 66cm | 50 - 55cm | 56cm | 60cm |
| XXS | S16FXXSM | S16SFCXXSM | S16CDLXXS | 66 - 76cm | 56 - 62cm | 64cm | 65cm |
| XS | S16FXSM | S16SFCXSM | S16CDLXS | 76 - 86cm | 60 - 66cm | 73cm | 70cm |
| S | S16FSM | S16SFCSM | S16CDLS | 86 - 96cm | 66 - 74cm | 80cm | 78cm |
| M | S16FMED | S16SFCMED | S16CDLM | 96 - 108cm | 74 - 80cm | 90cm | 86cm |
| L | S16FLRG | - | S16CDLL | 108 - 130cm | 80 - 88cm | 97cm | 94cm |
| XL | S16FXLRG | - | S16CDLXL | 130 - 144cm | 88 - 100cm | 114cm | 103cm |
| XXL | S16FXLRG | - | S16CDLXXL | 144 - 160cm | 96 - 104cm | 120cm | 107cm |

