

In-situ Straight Leg Sling

S17

- Designed to provide full body support on transfer
- Easy to apply with low-friction fabric to leg areas for ease of application
- Soft binding as standard
- Fixed colour coded loops or Silva-Safety Slot® tapes
- Once the patient has been transferred, this sling can be left in-situ, subject to an appropriate risk assessment
- Wash Temp 71°C/159°F
- Tumble dry low heat (max 50°C)



nb: Sling shown in image is with fixed Silva-Safety Slots® tapes

Product Codes & Specifications:



	Silva-Superfine PLUS®	Silva-Superfine® with COOLMAX®	Silva-Superfine® with COOLMAX® Double Layer	A	B	C	D
XS	-	S17SFCXSM	-	76 - 86cm	60 - 66cm	76cm	55cm
S	S17FSM	-	-	86 - 96cm	66 - 74cm	80cm	58cm
M	S17FMED	-	-	96 - 108cm	74 - 80cm	90cm	66cm
L	S17FLRG	-	-	108 - 130cm	80 - 88cm	102cm	74cm
XL	S17FXLRG	-	-	130 - 144cm	88 - 100cm	120cm	78cm

