

In-Situ High Easy

- Designed to provide full body & head support on transfer
- Adjustable integral head support with removable stiffeners
- Easy to apply with low-friction fabric to leg areas for ease of application
- Soft binding as standard
- Universal clip fixing
- Once the patient has been transferred, this sling can be left in-situ, subject to an appropriate risk assessment
- Wash Temp 71°C/159°F
- Tumble dry low heat (max 50°C)

HYAD



Product Codes & Specifications:



	Silva-Superfine PLUS®	Silva-Superfine® with COOLMAX®	Silva-Superfine® with COOLMAX® Double Layer	A	B	C	D
XXS	HYADFXXSMCP	HYADSFCXXSMCP	HYADSFCDLXXSMCP	66 - 76cm	56 - 62cm	66cm	64cm
XS	HYADFXSMCP	HYADSFXXSMCP	HYADSFCDLXSMP	76 - 86cm	60 - 66cm	68cm	69cm
S	HYADFSMCP	HYADSFCMCP	HYADSFCDLSMCP	86 - 96cm	66 - 74cm	76cm	77cm
M	HYADFMEDCP	HYADSFCMEDCP	HYADSFCDLMEDCP	96 - 108cm	74 - 80cm	86cm	84cm
L	HYADFLRGCP	-	HYADSFCDLLRGCP	108 - 130cm	80 - 88cm	98cm	93cm
XL	HYADFXLRGCP	-	HYADSFCDLXLRGCP	130 - 144cm	88 - 100cm	102cm	104cm
XXL	HYADFXLLRGCP	-	-	144 - 160cm	96 - 104cm	120cm	108cm

